

For Front or Rear Paw Knuckling

Walkin'® No-Knuckling Training Sock



A temporary training tool, for short-term, multiple use to help correct gait and improve paw placement. With specially designed socks for both front and rear paw knuckling. The No Knuckling Training Sock is designed to enhance proprioception in patients with neurological deficits and perfect for rehabilitative use for pets recovering from spinal surgery or suffering from disc disease. The sensory stimulation acts as a reminder for pets to lift their limb higher during gait training exercises.

Sizes: Front: XS – XL, Rear: XS – XL,


Retail Price: 45⁹⁵ – 69⁹⁵

Vet Professionals: Call for Special Pricing

Sizing for Front No Knuckling Training Sock

Size:	Top of paw to carpal:	Weight:
XSmall	2"	10 – 15 lbs.
Small	2.25"	15 – 25 lbs.
Medium	2.5"	25 – 60 lbs.
Large	2.75"	60 – 100 lbs.
X Large	3" or greater	Greater than 100 lbs.

Measure from top of paw to middle of carpal joint on your pet. Ideally try to do this with your pet in a standing position.




Advantages in a clinic setting:

- Utilize as a gait retraining tool
- Lightweight to target patients with limb weakness
- Designed to be used by rehab professionals as part of a treatment plan
- Instruct pet owners for use during at-home training sessions
- Use During: Enhanced Cavaletti Poles, Enhanced Figure 8's or Cone Weavings & Treadmill Patterning

Sizing for Rear No Knuckling Training Sock

Size:	Top of paw to hock:	Weight:
XSmall	Up to 3 inches	Up to 15 lbs.
Small	Up to 4 inches	Up to 30 lbs.
Medium	4 to 6 inches	30 – 60 lbs.
Large	6 to 8 inches	60 – 80 lbs.
X Large	8" or greater	Greater than 80 lbs.

Measure from top of paw to center of hock joint. Ideally try to do this with your pet in a standing position.



No-Knuckling Training Sock Reviews



No-Knuckling Training Sock is light and adjustable!

I was very surprised when I saw the No-Knuckling Training Sock in action today. As a neurologist we are always addressing two chronic concerns, dogs with knuckling postures and proprioceptive training. The No-Knuckling Training Sock is light and adjustable and really addresses both issues. Most of the other similar devices are heavier which presents a distinct disadvantage. Almost all of the other devices have to be custom-fitted but not the No-Knuckling Training Sock. The two sizes I reviewed will fit almost any size dog and that also will keep the cost down. I was very pleased with the device in use in a dog in our clinic. It represents a great improvement.

H. Steven Steinberg,
VMD, DACVIM (Neurology), CCRT
Chief of Staff VCA VRA
Instructor Canine Rehabilitation Institute,
Wellington, FL

Walking again after FCE!

I just used the No-Knuckling Training Sock on a 10-year-old lab that is relearning to walk after suffering from an FCE. He is very weak in both hind limbs and has trouble correctly placing his right hind foot, causing him to scrape his toes. I've tried other similar products, but because of their weight, it only made the dragging worse. The No-Knuckling Training Sock is very lightweight and a medium fit him perfectly. It was simple to put on and is easily adjustable. Even with an abnormal gait, the No-Knuckling Training Sock stayed in place and didn't move around. I will definitely use this product in the future.

April Milchak, CCRA

Great for Rehabilitation Patients

The No-Knuckling Training Sock has been very helpful for many of my rehabilitation patients. It is very lightweight and easy to put on. Also, I'm impressed with how comfortable it is. Dogs can wear it during numerous exercises or walks around the neighborhood and don't seem bothered by it at all.

Robin Price, RVT, CCRP